

For Research Use Only

Mouse Hif1a Recombinant antibody, PBS Only (Detector)



Catalog Number: 82989-2-PBS

Basic Information

Catalog Number:

82989-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33620

GenBank Accession Number:

NM_001313919

GeneID (NCBI):

15251

UNIPROT ID:

Q61221

Full Name:

hypoxia inducible factor 1, alpha
subunit

Calculated MW:

93 kDa

Purification Method:

Protein A purification

CloneNo.:

230101D11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00056-1, Mouse H1a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82989-3-PBS. Detection antibody: 82989-2-PBS. Standard: Ag33620. Range: 0.156-10 ng/mL.