

For Research Use Only

IL-1R1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82990-2-PBS



Basic Information

Catalog Number: 82990-2-PBS	GenBank Accession Number: BC075063	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3554	CloneNo.: 230192E12
Source: Rabbit	UNIPROT ID: P14778	
Isotype: IgG	Full Name: interleukin 1 receptor, type I	
Immunogen Catalog Number: AG25998	Calculated MW: 569 aa, 65 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00062-1, IL-1R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82990-2-PBS. Detection antibody: 82990-1-PBS. Standard: Ag25998. Range: 0.188-12 ng/mL