

For Research Use Only

# CD226 Recombinant antibody, PBS Only (Capture)

Catalog Number:82991-1-PBS



## Basic Information

Catalog Number:

82991-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074787

GeneID (NCBI):

10666

ENSEMBL Gene ID:

ENSG00000150637

UNIPROT ID:

Q15762

Full Name:

CD226 molecule

Calculated MW:

336 aa, 39 kDa

Purification Method:

Protein A purification

CloneNo.:

230041H8

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01106-1, Human CD226 Recombinant Matched Antibody Pair - PBS only. 82991-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard EG0027. 82991-2-PBS was HRP conjugated as the detection antibody. Range: 97.7-6250 pg/mL