

For Research Use Only

LAMP5 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82992-4-PBS

Basic Information

Catalog Number: 82992-4-PBS	GenBank Accession Number: BC050727	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 24141	CloneNo.: 230102H4
Source: Rabbit	UNIPROT ID: Q9UJQ1	
Isotype: IgG	Full Name: chromosome 20 open reading frame	
Immunogen Catalog Number: AG33894	103	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00069-1, LAMP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82992-4-PBS. Detection antibody: 82992-1-PBS. Standard: Ag33894. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00069-2, LAMP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82992-3-PBS. Detection antibody: 82992-4-PBS. Standard: Ag33894. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00069-3, LAMP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82992-4-PBS. Detection antibody: 82992-2-PBS. Standard: Ag33894. Range: 1.25-20 ng/mL