

For Research Use Only

Endoglin/CD105 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82993-2-PBS

Basic Information

Catalog Number:
82993-2-PBS

Size:
1 mg/ml

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC014271

GeneID (NCBI):
2022

UNIPROT ID:
P17813

Full Name:
endoglin

Calculated MW:
71 kDa

Purification Method:
Protein A purification

CloneNo.:
230047C5

Applications

Tested Applications:
IF, Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00089-1, ENG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82993-1-PBS. Detection antibody: 82993-2-PBS. Standard: Eg0587. Range: 0.156-10 ng/mL