

For Research Use Only

Endoglin/CD105 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82993-4-PBS

Basic Information

Catalog Number:

82993-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC014271

GeneID (NCBI):

2022

UNIPROT ID:

P17813

Full Name:

endoglin

Calculated MW:

71 kDa

Purification Method:

Protein A purification

CloneNo.:

230047H5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00089-2, ENG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82993-4-PBS. Detection antibody: 82993-3-PBS. Standard: Eg0587. Range: 0.156-10 ng/mL



Cytometric bead array standard curve of MP00089-2, ENG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82993-4-PBS. Detection antibody: 82993-3-PBS. Standard: Eg0587. Range: 2.5-20 ng/mL