

For Research Use Only

RILP Recombinant antibody, PBS Only (Capture)



Catalog Number: 82996-1-PBS

Basic Information

Catalog Number: 82996-1-PBS	GenBank Accession Number: BC031621	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 83547	CloneNo.: 230208C5
Source: Rabbit	UNIPROT ID: Q96NA2	Affinity: $K_D=1.30 \times 10^{-11}$
Isotype: IgG	Full Name: Rab interacting lysosomal protein	$K_{Off}=2.17 \times 10^{-5}$
Immunogen Catalog Number: AG4472	Calculated MW: 401 aa, 41 kDa	$K_{On}=1.67 \times 10^6$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00025-2, RLP Recombinant Matched Antibody Pair - PBS only. Capture antibody: 82996-1-PBS. Detection antibody: 82996-2-PBS. Standard: Ag4472. Range: 1.25-20 ng/mL