

For Research Use Only

FOXO1 Recombinant antibody, PBS Only (Detector)



Catalog Number:82997-3-PBS

Basic Information

Catalog Number: 82997-3-PBS	GenBank Accession Number: BC021981	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2308	CloneNo.: 230183H6
Source: Rabbit	UNIPROT ID: Q12778	
Isotype: IgG	Full Name: forkhead box O1	
Immunogen Catalog Number: AG13296	Calculated MW: 655 aa, 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00024-1, FOXO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82997-1-PBS. Detection antibody: 82997-3-PBS. Standard: Ag13296. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00024-2, FOXO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82997-2-PBS. Detection antibody: 82997-3-PBS. Standard: Ag13296. Range: 1.25-20 ng/mL