

For Research Use Only

# CCDC23 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82998-1-PBS



## Basic Information

Catalog Number:

82998-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33650

GenBank Accession Number:

BC029427

GeneID (NCBI):

374969

UNIPROT ID:

Q8N300

Full Name:

coiled-coil domain containing 23

Purification Method:

Protein A purification

CloneNo.:

230175B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00063-1, CCDC23 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82998-2-PBS. Detection antibody: 82998-1-PBS. Standard: Ag33650. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00063-2, CCDC23 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82998-3-PBS. Detection antibody: 82998-1-PBS. Standard: Ag33650. Range: 0.625-20 ng/mL.