

For Research Use Only

# UBN1-Specific Recombinant antibody, PBS Only (Detector)



Catalog Number:83001-1-PBS

## Basic Information

Catalog Number: 83001-1-PBS	GenBank Accession Number: NM_016936	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 29855	CloneNo.: 240014A8
Source: Rabbit	UNIPROT ID: Q9NPG3	
Isotype: IgG	Full Name: ubinuclain 1	
	Calculated MW: 121 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP00095-1, UBN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83001-2-PBS. Detection antibody: 83001-1-PBS. Standard: SY00809. Range: 0.078-5 ng/mL



Cytometric bead array standard curve of MP00095-2, UBN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83001-4-PBS. Detection antibody: 83001-1-PBS. Standard: SY00809. Range: 0.078-5 ng/mL



Cytometric bead array standard curve of MP00095-3, UBN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83001-3-PBS. Detection antibody: 83001-1-PBS. Standard: SY00809. Range: 0.078-5 ng/mL