

For Research Use Only

UBN1-Specific Recombinant antibody, PBS Only (Capture)



Catalog Number: 83001-4-PBS

Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 83001-4-PBS | GenBank Accession Number: NM_016936 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 29855 | CloneNo.: 240014H7 |
| Source: Rabbit | UNIPROT ID: Q9NPG3 | |
| Isotype: IgG | Full Name: ubinuclain 1 | |
| | Calculated MW: 121 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00095-2, UBN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83001-4-PBS. Detection antibody: 83001-1-PBS. Standard: SY00809. Range: 0.078-5 ng/mL.