

For Research Use Only

PIGM Recombinant antibody, PBS Only (Capture)



Catalog Number: 83006-1-PBS

Basic Information

Catalog Number: 83006-1-PBS	GenBank Accession Number: BC019865	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 93183	CloneNo.: 240025A6
Source: Rabbit	UNIPROT ID: Q9H3S5	
Isotype: IgG	Full Name: phosphatidylinositol glycan anchor biosynthesis, class M	
Immunogen Catalog Number: AG14245	Calculated MW: 423 aa, 49 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00048-1, PIGM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83006-1-PBS. Detection antibody: 83006-2-PBS. Standard: Ag14245. Range: 0.625-20 ng/mL