

For Research Use Only

PIGM Recombinant antibody, PBS Only (Detector)



Catalog Number: 83006-2-PBS

Basic Information

Catalog Number:

83006-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14245

GenBank Accession Number:

BC019865

GeneID (NCBI):

93183

UNIPROT ID:

Q9H3S5

Full Name:

phosphatidylinositol glycan anchor biosynthesis, class M

Calculated MW:

423 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

240025A10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00048-1, PIGM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83006-1-PBS. Detection antibody: 83006-2-PBS. Standard: Ag14245. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00048-2, PIGM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83006-3-PBS. Detection antibody: 83006-2-PBS. Standard: Ag14245. Range: 0.625-20 ng/mL