

For Research Use Only

IL12B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83008-2-PBS

Basic Information

Catalog Number:

83008-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC067498

GeneID (NCBI):

3593

UNIPROT ID:

P29460

Full Name:

interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40)

Calculated MW:

328 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

240008F6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00035-1, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-2-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00035-2, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-4-PBS. Detection antibody: 83008-2-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL