

For Research Use Only

# IL12B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83008-3-PBS

## Basic Information

<b>Catalog Number:</b> 83008-3-PBS	<b>GenBank Accession Number:</b> BC067498	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 3593	<b>CloneNo.:</b> 240008G6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P29460	
<b>Isotype:</b> IgG	<b>Full Name:</b> interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40)	
	<b>Calculated MW:</b> 328 aa, 37 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00035-1, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-2-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00035-3, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-1-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL