

For Research Use Only

golgin 97 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83009-1-PBS

Basic Information

Catalog Number: 83009-1-PBS	GenBank Accession Number: BC032853	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2800	CloneNo.: 230523A4
Source: Rabbit	UNIPROT ID: Q92805	
Isotype: IgG	Full Name: golgi autoantigen, golgin subfamily a, 1	
Immunogen Catalog Number: AG3335	Calculated MW: 767 aa, 88 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00031-1, GOLGA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83009-3-PBS. Detection antibody: 83009-1-PBS. Standard: Ag3335. Range: 0.625-20 ng/mL.