

For Research Use Only

# golgin 97 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83009-2-PBS

## Basic Information

Catalog Number:

83009-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3335

GenBank Accession Number:

BC032853

GeneID (NCBI):

2800

UNIPROT ID:

Q92805

Full Name:

golgi autoantigen, golgin subfamily  
a, 1

Calculated MW:

767 aa, 88 kDa

Purification Method:

Protein A purification

CloneNo.:

230523A6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00031-2, GOLGA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83009-3-PBS. Detection antibody: 83009-2-PBS. Standard: Ag3335. Range: 1.25-20 ng/mL.