

For Research Use Only

golgin 97 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83009-3-PBS

Basic Information

Catalog Number:

83009-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3335

GenBank Accession Number:

BC032853

GeneID (NCBI):

2800

UNIPROT ID:

Q92805

Full Name:

golgi autoantigen, golgin subfamily a, 1

Calculated MW:

767 aa, 88 kDa

Purification Method:

Protein A purification

CloneNo.:

230523G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00031-1, GOLGA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83009-3-PBS. Detection antibody: 83009-1-PBS. Standard: Ag3335. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00031-2, GOLGA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83009-3-PBS. Detection antibody: 83009-2-PBS. Standard: Ag3335. Range: 1.25-20 ng/mL.