

For Research Use Only

COL6A3-Specific Recombinant antibody, PBS Only (Capture)

Catalog Number: 83010-2-PBS



Basic Information

Catalog Number: 83010-2-PBS	GenBank Accession Number: NM_004369	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1293	CloneNo.: 230540C10
Source: Rabbit	UNIPROT ID: P12111	
Isotype: IgG	Full Name: collagen, type VI, alpha 3	
Immunogen Catalog Number: AG34752	Calculated MW: 344 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00027-2, COL6A3-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83010-2-PBS. Detection antibody: 83010-4-PBS. Standard: Ag34752. Range: 0.625-20 ng/mL.