

For Research Use Only

RLN2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83011-2-PBS

Basic Information

Catalog Number: 83011-2-PBS	GenBank Accession Number: BC126415	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6019	CloneNo.: 240058F11
Source: Rabbit	UNIPROT ID: P04090	
Isotype: IgG	Full Name: relaxin 2	
Immunogen Catalog Number: AG18781	Calculated MW: 185 aa, 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00040-1, RLN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83011-1-PBS. Detection antibody: 83011-2-PBS. Standard: Ag18781. Range: 1.25-20 ng/mL