

For Research Use Only

# GPR87 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83013-2-PBS

## Basic Information

Catalog Number:

83013-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23214

GenBank Accession Number:

BC132888

GeneID (NCBI):

53836

UNIPROT ID:

Q9BY21

Full Name:

G protein-coupled receptor 87

Purification Method:

Protein A purification

CloneNo.:

230453H2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00019-1, GPR87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83013-2-PBS. Detection antibody: 83013-3-PBS. Standard: Ag23214. Range: 0.625-20 ng/mL