

For Research Use Only

Mouse CD68 Recombinant antibody, PBS Only (Capture)



Catalog Number:83014-4-PBS

Basic Information

Catalog Number:

83014-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001291058.1

GeneID (NCBI):

12514

UNIPROT ID:

P31996-1

Full Name:

CD68 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230504H3

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00015-2 measured by cytometric bead array. Capture antibody: 83014-4-PBS. Detection antibody: 83014-2-PBS. Range: 0.625-20 ng/mL



Standard curve of MP00015-3 measured by cytometric bead array. Capture antibody: 83014-4-PBS. Detection antibody: 83014-1-PBS. Range: 0.625-20 ng/mL