

For Research Use Only

# PPP1R15B Recombinant antibody, PBS Only (Detector)

Catalog Number:83016-1-PBS



## Basic Information

<b>Catalog Number:</b> 83016-1-PBS	<b>GenBank Accession Number:</b> BC065280	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 84919	<b>CloneNo.:</b> 240052B7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q5SWA1	
<b>Isotype:</b> IgG	<b>Full Name:</b> protein phosphatase 1, regulatory (inhibitor) subunit 15B	
<b>Immunogen Catalog Number:</b> AG6213	<b>Calculated MW:</b> 79 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00002-1, PPP1R15B/CReP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83016-4-PBS. Detection antibody: 83016-1-PBS. Standard: Ag6213. Range: 0.625-20 ng/mL.