

For Research Use Only

PPP1R15B Recombinant antibody, PBS Only (Capture)

Catalog Number:83016-4-PBS



Basic Information

Catalog Number: 83016-4-PBS	GenBank Accession Number: BC065280	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84919	CloneNo.: 240052G6
Source: Rabbit	UNIPROT ID: Q5SWA1	
Isotype: IgG	Full Name: protein phosphatase 1, regulatory (inhibitor) subunit 15B	
Immunogen Catalog Number: AG6213	Calculated MW: 79 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00002-1, PPP1R15B/CReP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83016-4-PBS. Detection antibody: 83016-1-PBS. Standard: Ag6213. Range: 0.625-20 ng/mL.