

For Research Use Only

PPP1R15B Recombinant antibody, PBS Only (Capture)

Catalog Number:83016-5-PBS



Basic Information

Catalog Number:

83016-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6213

GenBank Accession Number:

BC065280

GeneID (NCBI):

84919

UNIPROT ID:

Q5SWA1

Full Name:

protein phosphatase 1, regulatory (inhibitor) subunit 15B

Calculated MW:

79 kDa

Purification Method:

Protein A purification

CloneNo.:

240052H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00002-3, PPP1R15B/CRP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83016-5-PBS. Detection antibody: 83016-2-PBS. Standard: Ag6213. Range: 0.625-20 ng/mL.