

For Research Use Only

ABCA3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83018-2-PBS

Basic Information

Catalog Number: 83018-2-PBS	GenBank Accession Number: BC140895	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 21	CloneNo.: 23019488
Source: Rabbit	UNIPROT ID: Q99758	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family A (ABC1), member 3	
Immunogen Catalog Number: AG24749		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00121-1, ABCA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83018-2-PBS. Detection antibody: 83018-3-PBS. Standard: Ag24749. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00121-2, ABCA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83018-2-PBS. Detection antibody: 83018-4-PBS. Standard: Ag24749. Range: 0.313-20 ng/mL