

For Research Use Only

ABCA3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83018-3-PBS

Basic Information

Catalog Number:

83018-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24749

GenBank Accession Number:

BC140895

GeneID (NCBI):

21

UNIPROT ID:

Q99758

Full Name:

ATP-binding cassette, sub-family A (ABC1), member 3

Purification Method:

Protein A purification

CloneNo.:

230194C1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00121-1, ABCA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83018-2-PBS. Detection antibody: 83018-3-PBS. Standard: Ag24749. Range: 0.313-20 ng/mL.