

For Research Use Only

# ABCA3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83018-5-PBS



## Basic Information

Catalog Number:

83018-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24749

GenBank Accession Number:

BC140895

GeneID (NCBI):

21

UNIPROT ID:

Q99758

Full Name:

ATP-binding cassette, sub-family A  
(ABC1), member 3

Purification Method:

Protein A purification

CloneNo.:

230194A9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00121-3, Human ABCA3 Recombinant Matched Antibody Pair - PBS only. 83018-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24749. 83018-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL

HeLa and SW480 cell lysates were measured. The human ABCA3 concentration of detected samples was determined to be 2.8 ng/mL (based on a 2.1 mg/mL extract load) in HeLa cell lysate and 1.9 ng/mL (based on a 2.3 mg/mL extract load) in SW480 cell lysate.