

For Research Use Only

Cadherin-19 Recombinant antibody, PBS Only (Capture)

Catalog Number:83019-4-PBS



Basic Information

Catalog Number: 83019-4-PBS	GenBank Accession Number: BC015877	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 28513	CloneNo.: 230271F2
Source: Rabbit	UNIPROT ID: Q9H159	
Isotype: IgG	Full Name: cadherin 19, type 2	
Immunogen Catalog Number: AG2000	Calculated MW: 55 kDa; 87 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00044-1, Cadherin-19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83019-4-PBS. Detection antibody: 83019-3-PBS. Standard: Ag2000. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00044-3, Cadherin-19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83019-4-PBS. Detection antibody: 83019-2-PBS. Standard: Ag2000. Range: 0.625-20 ng/mL