

For Research Use Only

ANO9 Recombinant antibody, PBS Only (Detector)



Catalog Number:83025-1-PBS

Basic Information

Catalog Number: 83025-1-PBS	GenBank Accession Number: NM_001012302	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 338440	CloneNo.: 230256A7
Source: Rabbit	UNIPROT ID: A1A5B4	
Isotype: IgG	Full Name: anoctamin 9	
	Calculated MW: 90 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00092-1, ANO9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83025-5-PBS. Detection antibody: 83025-1-PBS. Standard: SY00332. Range: 0.078-5 ng/mL