

For Research Use Only

# ANO9 Recombinant antibody, PBS Only (Capture)



Catalog Number:83025-6-PBS

## Basic Information

Catalog Number:

83025-6-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001012302

GeneID (NCBI):

338440

UNIPROT ID:

A1A5B4

Full Name:

anoctamin 9

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

230256E10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00092-3, ANO9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83025-6-PBS. Detection antibody: 83025-2-PBS. Standard: SY00332. Range: 0.078-5 ng/mL