

For Research Use Only

# ACVR1B Recombinant antibody, PBS Only (Capture)



Catalog Number:83026-1-PBS

## Basic Information

<b>Catalog Number:</b> 83026-1-PBS	<b>GenBank Accession Number:</b> BC000254	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 91	<b>CloneNo.:</b> 230260C5
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P36896	
<b>Isotype:</b> IgG	<b>Full Name:</b> activin A receptor, type IB	
<b>Immunogen Catalog Number:</b> AG0129	<b>Calculated MW:</b> 57 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00098-2, ACVR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83026-1-PBS. Detection antibody: 83026-4-PBS. Standard: Ag0129. Range: 0.313-20 ng/mL.