

For Research Use Only

ACVR1B Recombinant antibody

Catalog Number: 83026-6-RR



Basic Information

Catalog Number: 83026-6-RR	GenBank Accession Number: BC000254	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 91	CloneNo.: 230260C1
Source: Rabbit	UNIPROT ID: P36896	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG	Full Name: activin A receptor, type 1B	
Immunogen Catalog Number: AG0129	Calculated MW: 57 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Caco-2 cells, HepG2 cells, HEK-293 cells
Species Specificity: Human	

Background Information

ACVR1B (ALK4, activin A receptor type 1B) is a transforming growth factor-beta (TGF-beta) superfamily member (PMID: 37020544). A potential role for ACVR1B in regulating muscle mass is that it is part of the transforming growth factor b (TGFb) pathway regulating myostatin (PMID: 21063444).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83026-6-RR (ACVR1B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.