

For Research Use Only

ACVR1B Recombinant antibody, PBS Only (Detector)

Catalog Number:83026-7-PBS



Basic Information

Catalog Number: 83026-7-PBS	GenBank Accession Number: BC000254	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 91	CloneNo.: 230260D10
Source: Rabbit	UNIPROT ID: P36896	
Isotype: IgG	Full Name: activin A receptor, type IB	
Immunogen Catalog Number: AG0129	Calculated MW: 57 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00098-4, Human ACVR1B Recombinant Matched Antibody Pair - PBS only. 83026-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0129. 83026-7-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL



The mean ACVR1B concentration was determined to be 40.70 ng/mL in HepG2 cell extract based on a 1.00 mg/mL extract load.