

For Research Use Only

CLDN18 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83028-1-PBS

Basic Information

Catalog Number: 83028-1-PBS	GenBank Accession Number: NM_001002026	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51208	CloneNo.: 230198A10
Source: Rabbit	UNIPROT ID: P56856-2	
Isotype: IgG	Full Name: claudin 18	
Immunogen Catalog Number: AG34049	Calculated MW: 28 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00105-1, CLDN18 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83028-2-PBS. Detection antibody: 83028-1-PBS. Standard: Ag34049. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00105-2, CLDN18 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83028-1-PBS. Detection antibody: 83028-3-PBS. Standard: Ag34049. Range: 0.313-20 ng/mL.