

For Research Use Only

CLDN18 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83028-2-PBS

Basic Information

Catalog Number:	83028-2-PBS	GenBank Accession Number:	NM_001002026	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	51208	CloneNo.:	230198B4
Source:	Rabbit	UNIPROT ID:	P56856-2		
Isotype:	IgG	Full Name:	claudin 18		
Immunogen Catalog Number:	AG34049	Calculated MW:	28 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00105-1, CLDN18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83028-2-PBS. Detection antibody: 83028-1-PBS. Standard: Ag34049. Range: 0.313-20 ng/mL.