

For Research Use Only

BMPR1B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83029-1-PBS

Basic Information

Catalog Number: 83029-1-PBS	GenBank Accession Number: BC047773	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 658	CloneNo.: 230244B11
Source: Rabbit	UNIPROT ID: O00238	Affinity: $K_D = 1.49 \times 10^{-10}$
Isotype: IgG	Full Name: bone morphogenetic protein receptor, type IB	$K_{Off} = 1.22 \times 10^{-4}$ $K_{On} = 8.14 \times 10^4$
Immunogen Catalog Number: AG24902	Calculated MW: 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00108-2, BMPRII Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83029-3-PBS. Detection antibody: 83029-1-PBS. Standard: Ag24902. Range: 0.313-10 ng/mL.