

For Research Use Only

BMPR1B Recombinant antibody, PBS Only (Detector)



Catalog Number: 83029-2-PBS

Basic Information

| | | |
|---|---|--|
| Catalog Number: 83029-2-PBS | GenBank Accession Number: BC047773 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 658 | CloneNo.: 230244D8 |
| Source: Rabbit | UNIPROT ID: O00238 | Affinity: $K_D = 4.61 \times 10^{-10}$ |
| Isotype: IgG | Full Name: bone morphogenetic protein receptor, type IB | $K_{Off} = 2.71 \times 10^{-4}$ $K_{On} = 5.89 \times 10^4$ |
| Immunogen Catalog Number: AG24902 | Calculated MW: 57 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00108-1, BMPRII Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83029-7-PBS. Detection antibody: 83029-2-PBS. Standard: Ag24902. Range: 0.313-10 ng/mL