

For Research Use Only

BMPR1B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83029-5-PBS

Basic Information

Catalog Number:	83029-5-PBS	GenBank Accession Number:	BC047773	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	658	CloneNo.:	230244F3
Source:	Rabbit	UNIPROT ID:	O00238		
Isotype:	IgG	Full Name:	bone morphogenetic protein receptor, type IB		
Immunogen Catalog Number:	AG24902	Calculated MW:	57 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00108-3, BMPRII Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83029-5-PBS. Detection antibody: 83029-4-PBS. Standard: Ag24902. Range: 0.625-10 ng/mL