

For Research Use Only

# GABRA4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83030-1-PBS

## Basic Information

Catalog Number:

83030-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3600

GenBank Accession Number:

BC035055

GeneID (NCBI):

2557

UNIPROT ID:

P48169

Full Name:

gamma-aminobutyric acid (GABA) A receptor, alpha 4

Calculated MW:

554 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

230302C6

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00115-1 measured by cytometric bead array. Capture antibody: 83030-2-PBS. Detection antibody: 83030-1-PBS. Range: 0.313-20 ng/mL



Standard curve of MP00115-3 measured by cytometric bead array. Capture antibody: 83030-3-PBS. Detection antibody: 83030-1-PBS. Range: 1.25-20 ng/mL