

For Research Use Only

SLC29A3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83036-1-PBS

Basic Information

Catalog Number:

83036-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18563

GenBank Accession Number:

BC120996

GeneID (NCBI):

55315

UNIPROT ID:

Q9BZD2

Full Name:

solute carrier family 29 (nucleoside transporters), member 3

Calculated MW:

475 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

230281C12

Applications

Tested Applications:

WB, Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with 83036-1-RR (SLC29A3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83036-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00061-1, SLC29A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83036-2-PBS. Detection antibody: 83036-1-PBS. Standard: Ag18563. Range: 0.625-20 ng/mL



Bi-layer interferometry (BLI) kinetic assays of 83036-1-RR against Human SLC29A3 were performed. The affinity constant is 0.455 nM.