

For Research Use Only

# SLC29A3 Recombinant antibody, PBS Only (Capture)



Catalog Number:83036-2-PBS

## Basic Information

|   |   |   |
|---|---|---|
| <b>Catalog Number:</b><br>83036-2-PBS       | <b>GenBank Accession Number:</b><br>BC120996                                      | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1mg/ml                      | <b>GeneID (NCBI):</b><br>55315  | <b>CloneNo.:</b><br>230281E12                         |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q9BZD2  |   |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>solute carrier family 29 (nucleoside transporters), member 3 |   |
| <b>Immunogen Catalog Number:</b><br>AG18563 | <b>Calculated MW:</b><br>475 aa, 52 kDa   |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00061-1, SLC29A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83036-2-PBS. Detection antibody: 83036-1-PBS. Standard: Ag18563. Range: 0.625-20 ng/mL