

For Research Use Only

DRD1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83037-5-PBS



Basic Information

Catalog Number: 83037-5-PBS	GenBank Accession Number: BC074978	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1812	CloneNo.: 241782B12
Source: Rabbit	UNIPROT ID: P21728	
Isotype: IgG	Full Name: dopamine receptor D1	
Immunogen Catalog Number: AG12366	Calculated MW: 446 aa, 49 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01326-1, DRD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83037-5-PBS. Detection antibody: 83037-7-PBS. Standard: Ag12366. Range: 0.313-40 ng/mL