

For Research Use Only

ABCD1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83038-2-PBS

Basic Information

Catalog Number:

83038-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0453

GenBank Accession Number:

BC025358

GeneID (NCBI):

215

UNIPROT ID:

P33897

Full Name:

ATP-binding cassette, sub-family D (ALD), member 1

Calculated MW:

745 aa, 83 kDa

Purification Method:

Protein A purification

CloneNo.:

230226F6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 78-5000 pg/mL.



Cytometric bead array standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 0.625-20 ng/mL.