

For Research Use Only

ABCD1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83038-3-PBS

Basic Information

Catalog Number: 83038-3-PBS	GenBank Accession Number: BC025358	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 215	CloneNo.: 230226G11
Source: Rabbit	UNIPROT ID: P33897	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family D (ALD), member 1	
Immunogen Catalog Number: AG0453	Calculated MW: 745 aa, 83 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP00070-1, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-1-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 78-5000 pg/mL



Sandwich ELISA standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 78-5000 pg/mL



Cytometric bead array standard curve of MP00070-1, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-1-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00070-2, ABCD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83038-2-PBS. Detection antibody: 83038-3-PBS. Standard: Ag0453. Range: 0.625-20 ng/mL