

For Research Use Only



Acetyl-Histone H2A (Lys9) Recombinant antibody, PBS Only (Detector)

Catalog Number: **83041-2-PBS**

Basic Information

Catalog Number: 83041-2-PBS	GenBank Accession Number: BC093836	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3012	CloneNo.: 230349C9
Source: Rabbit	UNIPROT ID: P04908	
Isotype: IgG	Full Name: histone cluster 1, H2ae	
	Calculated MW: 14 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00085-3, HIST1H2AE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83041-1-PBS. Detection antibody: 83041-2-PBS. Standard: SY02140Ac. Range: 0.156-5 ng/mL