

For Research Use Only

# HIST1H2AE Recombinant antibody, PBS Only (Detector)



Catalog Number:83041-3-PBS

## Basic Information

Catalog Number:

83041-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093836

GeneID (NCBI):

3012

UNIPROT ID:

P04908

Full Name:

histone cluster 1, H2ae

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

230349D4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00085-2, HIST1H2AE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83041-1-PBS. Detection antibody: 83041-3-PBS. Standard: SY02140Ac. Range: 0.156-5 ng/mL