

For Research Use Only

HIST1H2AE Recombinant antibody, PBS Only (Detector)



Catalog Number:83041-4-PBS

Basic Information

Catalog Number:

83041-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093836

GeneID (NCBI):

3012

UNIPROT ID:

P04908

Full Name:

histone cluster 1, H2ae

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

230349E5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00085-1, HIST1H2AE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83041-1-PBS. Detection antibody: 83041-4-PBS. Standard: SY02140Ac. Range: 0.156-5 ng/mL