

For Research Use Only

MTA3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83042-2-PBS

Basic Information

Catalog Number: 83042-2-PBS	GenBank Accession Number: BC053631	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 57504	CloneNo.: 230350B12
Source: Rabbit	UNIPROT ID: Q9BTC8	
Isotype: IgG	Full Name: metastasis associated 1 family, member 3	
Immunogen Catalog Number: AG6400	Calculated MW: 68 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00017-1, MTA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83042-4-PBS. Detection antibody: 83042-2-PBS. Standard: Ag6400. Range: 0.625-20 ng/mL



Biolayer interferometry (BLI) kinetic assays of 83042-2-RR against Human MTA3 were performed. The affinity constant is below 1 pM.