

For Research Use Only

# TLR3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83043-2-PBS

## Basic Information

<b>Catalog Number:</b> 83043-2-PBS	<b>GenBank Accession Number:</b> BC094737	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 7098	<b>CloneNo.:</b> 230253H6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> O15455	
<b>Isotype:</b> IgG	<b>Full Name:</b> toll-like receptor 3	
<b>Immunogen Catalog Number:</b> AG12171	<b>Calculated MW:</b> 904 aa, 104 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00114-1 measured by cytometric bead array. Capture antibody: 83043-2-PBS. Detection antibody: 83043-1-PBS. Range: 0.625-20 ng/mL